



Meer info? Plus d'info? More info?

T: +32 3 870 60 80

E: advice@gillain.com

www.gillain.com

Size comparison of KEOFITT Valves

Valve	Typical viscosity	Particle diameter	Typical products
4KAI	1 – 100 cP	< 1,5 mm	Water, soft-drinks
W9	1 – 1.000 cP	< 3 mm	Beer, juice, milk products, pharm. products
152	1 – 800 cP	< 2 mm	Cells, blood products, insulin
W15	1 – 50.000 cP	< 8 mm	Yogurt, ketchup, honey, creams
W25	1 – 250.000 cP	< 17 mm	Fruit pulp, any liquids containing large particles
Basix	1 – 100 cP	< 1,5 mm	Water
Sesame	1 – 1.000 cP	< 3 mm	Beer, drinks, milk products, pharm. products
Micro Port	1 – 50 cP	< 0,5 mm	Filtered beer, milk, water, pharm. products

Membranes







icone

...-



